



REF CDC-05531-PK1A (IPNXXXXXX)

3 Fr. 1 Lumen 55 cm catheter length .016 inch dia. spring-wire guide

One-Lumen PICC



Contents:

- 1: One Lumen Catheter: 3 Fr. (1.03 mm OD) x 55 cm, T-Port Connector, Contamination Guard and Placement Wire Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: GlideThru™ Peel-Away Sheath: 3 Fr. x 2-3/4" (7 cm) Radiopaque over 3 Fr. Dilator Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Spring Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip On Both Ends) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Syringe: 10 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
2: 10 mL Luer-Lock Pre-Filled Saline Syringe Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®2 One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
- 1: Dust Cap: Non-Vented Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: SharpsAway® Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Catheter Trimmer Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Towel: 24" x 36" (61 cm x 91 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Towel Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Sterile Procedure Sign Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Patient ID Card Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Chart Sticker Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Patient Information Booklet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Checklist/CLIP Sheet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
2: Paper Tape Measure Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Tourniquet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX

- 5: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
10: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Surgical Apparel: Impervious Gown Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Dressing: BIOPATCH®3 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Surgical Apparel: Bouffant Cap Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Dressing: STATLOCK®4 PICC Plus Catheter Stabilization Device Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Surgical Apparel: Mask Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Surgical Apparel: Mask with Eye Shield Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: Tape: Steri-Strip®5 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
1: HemoHopper® Fluid Receptacle Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX
- 1A trademark of Becton, Dickinson and Company.
2A registered trademark of CareFusion or one of its subsidiaries.
3A registered trademark of Johnson & Johnson Corporation.
4A registered trademark of C. R. Bard, Inc.
5A registered trademark of 3M Company.

Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.
Contraindications: The PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)
Distal (20 Ga.)	0.3	103	1044

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.
†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

