



REF CA-04220 (IPNXXXXXX)

20^{Ga.} 3.81^{cm} catheter length .018^{inch dia. spring-wire guide}

QuickFlash® Radial Artery Catheterization

Contents:

- 1: QuickFlash® Radial Artery Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
 - 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm)
 - 1: Syringe: 3 mL Luer-Lock
 - 1: SharpsAway® Disposal Cup
 - 1: SharpsAway® II Locking Disposal Cup
 - 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
 - 1: Towel: 13" x 18" (33 cm x 46 cm)
 - 1: Needle Holder
 - 1: Safety Scalpel: #11
 - 1: Medication Label: 1% Lidocaine
 - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Dressing: I.V. Advanced 10 cm x 12 cm
 - 1: Surgical Apparel: Bouffant Cap
 - 1: Surgical Apparel: Mask with Eye Shield
 - 1: Suture: 3-0 Silk with Curved Needle
 - 1: Tubing±: 6" Extension with One-Way Stopcock
- ¹A trademark of Becton, Dickinson and Company.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
±Precaution: Contains Phthalate: DEHP.
Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.
Not made with natural rubber latex.
Fluid path components are non-pyrogenic.
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA



Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com
LBL080463 R00 (2025-09)

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



REF CA-04220

20^{Ga.} 3.81^{cm} catheter length .018^{inch dia. spring-wire guide}



QuickFlash® Radial Artery Catheterization