

REF **AK-09000** (IPNXXXXXX)

Peritoneal Lavage

8^{Fr.}

24 cm catheter length

.035 inch dia.
spring-wire
guide

Contents:

- 1: Peritoneal Lavage Catheter: 8 Fr. x 9-1/4" (24 cm) Radiopaque FEP with Side Holes Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Injection Needle: Eclipse™ 1 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Injection Needle: SafetyGlide™ 1 23 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 2: Syringe: 5 mL Luer-Slip Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Tubing Sets: Fluid Administration Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

¹A trademark of Becton, Dickinson and Company.

²A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm

⚠ WARNING: Cancer and
www.P65Warnings.ca.gov

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction.

Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in

the case of long-term exposure or application. Medical products contain

temporarily with pregnant women, i

Not made with natural rubber latex.

Store between 20 – 25°C. (68 – 77°F)

Sterilized using ethylene oxide unless otherwise specified.

Fluid path components are non-pyrogenic

Arrow International provides the enclosed medication label(s).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



LOT Lot Number



(01)XXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Set Product of xxx
Packaged in xxx

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL080525 R00 (2025-09)

REF AK-09000

8^{Fr.}

24 cm catheter length

.035 inch dia.
spring-wire
guide



Peritoneal Lavage

