



CAN-04001-AX

(IPNXXXXXX)



ARROW®

.018<sup>inch dia.
spring-wire
guide</sup>

Access Tray with TipTracker™ Stylet for use with 5 Fr. PICC



Contents:

- 1: GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator

1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer

1: VPS TipTracker™ Stylet with T-Port Connector: 0.015" (0.38 mm) dia. x 29" (74 cm)

1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm)

1: Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip

1: Syringe: 3 mL Luer-Lock

1: Dust Cap: Non-Vented

1: SharpsAway® II Locking Disposal Cup

1: SharpsAway® Disposal Cup

1: Catheter Trimmer

1: Towel: 24" x 36" (61 cm x 91 cm)

1: Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations

2: Towel

1: Needle Holder

1: Sterile Procedure Sign

2: Paper Tape Measure

2: Tourniquet

1: Foam Electrodes (3 per pack)

5: Gauze Pad: 2" x 2" (5 cm x 5 cm)

10: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Surgical Apparel: Impervious Gown

1: Dressing: I.V. Advanced 10 cm x 15.5 cm

1: Transducer Cover with Gel and Rubber Bands

1: Remote Cover

1: Surgical Tape: Transpore™² 1" x 54" Roll

1: Valve: MaxZero™³ Needleless Connector

1: T-Piece Cover

¹A trademark of Becton, Dickinson and Company.

²A trademark of 3M Company.

³A trademark of CareFusion or one of its subsidiaries.

Not made with natural rubber latex.
Store below 30°C (86°F). Avoid moisture, direct heat and direct light.

