



REF CDC-35563-HPK1A (IPNXXXXXX)

6 Fr. 3 Lumen 55 cm catheter length .018 inch dia. spring-wire guide

ARROW

ErgoPack

UDI



Pressure Injectable Three-Lumen PICC

Contents:

- | | | |
|---|--|---|
| 1: Three Lumen Catheter: 6 Fr. (2.14 mm OD) x 55 cm, Pressure Injectable, T-Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXX | 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX |
| 1: GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 3: Dust Cap: Non-Vented Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Checklist/CLIP Sheet Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX |
| 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: SecondSite™ Adjustable Hub: Fastener Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 2: Paper Tape Measure Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX |
| 1: Injection Needle: Eclipse™ 1 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: SecondSite™ Adjustable Hub: Catheter Clamp Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Tourniquet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 1: Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 5: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: SharpsAway® Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 10: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Catheter Trimmer Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Surgical Apparel: Impervious Gown Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 1: Syringe: 10 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Towel: 24" x 36" (61 cm x 91 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Dressing: BIOPATCH® Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 2: 10 mL Luer-Lock Pre-Filled Saline Syringe Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Towel Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Dressing: STATLOCK® PICC Plus Catheter Stabilization Device Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Surgical Apparel: Bouffant Cap Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Sterile Procedure Sign Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Surgical Apparel: Mask Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| | 1: Patient ID Card Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Surgical Apparel: Mask with Eye Shield Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| | 1: Chart Sticker Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Tape: Steri-Strip® Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |
| | 1: Patient Information Booklet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: HemoHopper® Fluid Receptacle Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX |

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

Pump
Flow Rate††
(mL/hr)

MAX Pressure
Injection
Flow Rate**
(mL/sec)

Distal (17 Ga.)	0.68	938	5950	6
Medial (18 Ga.)	0.49	277	2540	not rated
Proximal (18 Ga.)	0.5	284	2500	not rated

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.

