



ARROW®

REF AK-04020 (IPNXXXXXX)

20^{Ga.} 4.45^{cm} catheter length .018^{inch dia.} spring-wire guide

Radial Artery Catheterization


- Contents:**
- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
 - 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
 - 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
 - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
 - 1: Syringe: 3 mL Luer-Lock
 - 1: Foil Pack 70% Alcohol Prep Pad
 - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
 - 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
 - 1: Wing Clip: Suture
 - 1: Towel: 13" x 18" (33 cm x 46 cm)
 - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL082294 R00 (2026-05)



ARROW®

REF AK-04020 (IPNXXXXXX)

20^{Ga.} 4.45^{cm} catheter length .018^{inch dia.} spring-wire guide

Radial Artery Catheterization


- Contents:**
- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
 - 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
 - 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
 - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
 - 1: Syringe: 3 mL Luer-Lock
 - 1: Foil Pack 70% Alcohol Prep Pad
 - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
 - 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
 - 1: Wing Clip: Suture
 - 1: Towel: 13" x 18" (33 cm x 46 cm)
 - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL082294 R00 (2026-05)



REF AK-04020

ARROW®

20^{Ga.} 4.45^{cm} catheter length .018^{inch dia.} spring-wire guide

Radial Artery Catheterization



REF AK-04020

ARROW®

20^{Ga.} 4.45^{cm} catheter length .018^{inch dia.} spring-wire guide

Radial Artery Catheterization