



ARROW®

REF

AK-04020-C

(IPNXXXXXX)

20^{Ga.}

4.45^{cm}
catheter
length

.018^{inch dia.}
spring-wire
guides

Radial Artery Catheterization

- Contents:
- 1:

Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1:

Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1:

Injection Needle: 27 Ga. x 1/2" (1.27 cm)
- 1:

Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1:

Syringe: 3 mL Luer-Lock
- 1:

Foil Pack 70% Alcohol Prep Pad
- 1:

5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1:

1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
- 1:

Wing Clip: Suture
- 1:

Towel: 13" x 18" (33 cm x 46 cm)
- 1:

Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1:

Scalpel: #11
- 4:

Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4:


Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1:

Suture: 3-0 Braided Silk with Straight Needle
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.



LOT

Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX

(17)YYMMDD

(11)YYMMDD

(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx

Packaged in xxx

Teleflex®

www.teleflex.com

LBL082295 R00 (2026-05)



ARROW®

REF

AK-04020-C

(IPNXXXXXX)

20^{Ga.}

4.45^{cm}
catheter
length

.018^{inch dia.}
spring-wire
guides

Radial Artery Catheterization

- Contents:
- 1:

Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1:

Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1:

Injection Needle: 27 Ga. x 1/2" (1.27 cm)
- 1:

Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1:

Syringe: 3 mL Luer-Lock
- 1:

Foil Pack 70% Alcohol Prep Pad
- 1:

5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1:

1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
- 1:

Wing Clip: Suture
- 1:

Towel: 13" x 18" (33 cm x 46 cm)
- 1:

Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1:

Scalpel: #11
- 4:

Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4:


Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1:

Suture: 3-0 Braided Silk with Straight Needle
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.



LOT

Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX

(17)YYMMDD

(11)YYMMDD

(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx

Packaged in xxx

Teleflex®

www.teleflex.com

LBL082295 R00 (2026-05)



REF

AK-04020-C

ARROW®

20^{Ga.}

4.45^{cm}
catheter
length

.018^{inch dia.}
spring-wire
guides

Radial Artery Catheterization



REF

AK-04020-C

ARROW®

20^{Ga.}

4.45^{cm}
catheter
length

.018^{inch dia.}
spring-wire
guides

Radial Artery Catheterization