



## ARROWgard Blue<sup>®</sup> Percutaneous Sheath Introducer Kit Cath-Gard<sup>®</sup> for use with 7.5 - 8 Fr. Catheters

REF CDC-29903-X1A

9 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide

ErgoPack<sup>®</sup>

**Contents:**

- 1: Sheath: 9 Fr. x 10 cm Arrow-Flex<sup>®</sup> Radiopaque Polyurethane with ARROWgard<sup>®</sup> Antimicrobial Surface Treatment<sup>1</sup>, Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow<sup>®</sup> Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide<sup>™2</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: Injection Needle: SafetyGlide<sup>™2</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: 5 mL 1% Lidocaine HCl Solution
- 2: 10 mL Luer-Lock Pre-Filled Saline Syringe

- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep<sup>®3</sup> One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint (non-sterile solution)
- 1: Obturator: 8 Fr.
- 1: Cath-Gard<sup>®</sup> Catheter Contamination Shield: 80 cm with TwistLock<sup>™</sup> Adapter
- 1: Dust Cap: Non-Vented
- 1: SharpsAway<sup>®</sup> Disposal Cup
- 1: SharpsAway II<sup>™</sup> Locking Disposal Cup
- 1: Maximal Barrier Drape<sup>™</sup> with 4" fenestration
- 1: Towel
- 1: Local Anesthetic Component
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign

- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm<sup>®4</sup> 10 cm x 12 cm
- 1: Surgical Apparel: Impervious Gown
- 1: Dressing: BIOPATCH<sup>®5</sup>
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tubing<sup>3</sup>: 20" (50 cm) Extension
- 1: HemoHopper<sup>®</sup> Fluid Receptacle

<sup>1</sup>Licensed under US Patent No. 6,706,024.  
<sup>2</sup>A trademark of Becton, Dickinson and Company.  
<sup>3</sup>A registered trademark of CareFusion or one of its subsidiaries.  
<sup>4</sup>A registered trademark of 3M Company.  
<sup>5</sup>A registered trademark of Johnson & Johnson Corporation.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

**Contraindications: The ARROWgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

~~Latex~~  
 Not made with natural rubber latex

PHT DEHP  
 †Precaution: Contains Phthalate: DEHP.

20°C (68°F) - 25°C (77°F)  
 Store between 20 - 25°C (68 - 77°F).