



ARROWgard Blue[®] PSI Kit Cath-Gard[®] for use with 7.5 - 8 Fr. Catheters

ARROW[®]

REF **CDC-29903-1A**

9 Fr. | **10** cm sheath length | **.035** inch dia. spring-wire guide

ErgoPack[®]

Contents:

- 1: Sheath: 9 Fr. x 10 cm Arrow-Flex[®] Radiopaque Polyurethane with ARROWgard[®] Antimicrobial Surface Treatment¹, Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide^{™2} 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide^{™2} 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution

- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep^{®3} One-Step Solution with Hi-Lite Orange[™] Tint (non-sterile solution)
- 1: Obturator: 8 Fr.
- 1: Cath-Gard[®] Catheter Contamination Shield: 80 cm with TwistLock[™] Adapter
- 1: Dust Cap: Non-Vented
- 1: SharpsAway[®] Disposal Cup
- 1: SharpsAway II[™] Locking Disposal Cup
- 1: Maximal Barrier Drape[™] with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)

- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tubing[‡]: 20" (50 cm) Extension
- 1: HemoHopper[®] Fluid Receptacle

¹Licensed under US Patent No. 6,706,024.
²A trademark of Becton, Dickinson and Company.
³A registered trademark of CareFusion or one of its subsidiaries.



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.
Fluid path components are non-pyrogenic.

Contraindications: The ARROWgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Manufacturer:
Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex[®]
www.teleflex.com

S-29903-201D, Rev. 01 (2021-01)

DuPont™ Tyvek[®] packaging



ARROWgard Blue[®] PSI Kit Cath-Gard[®] for use with 7.5 - 8 Fr. Catheters

ARROW[®]

REF **CDC-29903-1A**

ErgoPack[®]

9 Fr. | **10** cm sheath length | **.035** inch dia. spring-wire guide