

Dimensions:
Width: 4" (10.2 cm)
Length: 11" (27.9 cm)

DO NOT BEND

REF **SAC-00324**

Arterial
Catheterization Set

ARROW

- cs** Arteriální katetrizační set
- da** Arterielt kateteriseringssæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:

- 1: Indwelling Catheter: 24 Ga x 2.5 cm Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm)
(Straight Soft Tip on Both Ends)
- 1: Catheter: 24 Ga x 3/4" (1.9 cm) Radiopaque over 26 Ga Introducer Needle

Catheter Priming Volume: 0.14 ml

Priming Volume is approximate.

Contraindications: Inadequate collateral arterial circulation.

Fluid path components are non-pyrogenic.



Arrow International, Inc., Subsidiary of Teleflex Incorporated
2400 Bernville Road, Reading, PA 19605 USA

EC REP Teleflex Medical, IDA Business and Technology Park
Dublin Road, Athlone, Co. Westmeath, Ireland

Teleflex
www.teleflex.com

LOT



SZ-00324-101D (3/17)

ARROW



24 Ga **2.5** cm catheter length

DO NOT BEND

DO NOT BEND

REF **SAC-00324**

Arterial
Catheterization Set

ARROW

- cs** Arteriální katetrizační set
- da** Arterielt kateteriseringssæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:

- 1: Indwelling Catheter: 24 Ga x 2.5 cm Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm)
(Straight Soft Tip on Both Ends)
- 1: Catheter: 24 Ga x 3/4" (1.9 cm) Radiopaque over 26 Ga Introducer Needle

Catheter Priming Volume: 0.14 ml

Priming Volume is approximate.

Contraindications: Inadequate collateral arterial circulation.

Fluid path components are non-pyrogenic.



Arrow International, Inc., Subsidiary of Teleflex Incorporated
2400 Bernville Road, Reading, PA 19605 USA

EC REP Teleflex Medical, IDA Business and Technology Park
Dublin Road, Athlone, Co. Westmeath, Ireland

Teleflex
www.teleflex.com

LOT



SZ-00324-101D (3/17)

ARROW



24 Ga **2.5** cm catheter length

DO NOT BEND