

DO NOT BEND

ARROW



REF SAC-00522-PBX

22 Ga. | 5 cm catheter length

Arterial Catheterization Set

- cs** Arteriální katetrizační set
- da** Arterielt kateteriseringssæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:

1: Indwelling Catheter: 22 Ga. x 5 cm Radiopaque PEBA

1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm)
(Straight Soft Tip on Both Ends)

1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW

Catheter Priming Volume: 0.16 mL

Priming Volume is approximate.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

⚠ California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

SZ-00522-103B (3/18)

LOT



REF SAC-00522-PBX

22 Ga. | 5 cm catheter length

Arterial Catheterization Set

DO NOT BEND

ARROW

DO NOT BEND

ARROW



REF SAC-00522-PBX

22 Ga. | 5 cm catheter length

Arterial Catheterization Set

- cs** Arteriální katetrizační set
- da** Arterielt kateteriseringssæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:

1: Indwelling Catheter: 22 Ga. x 5 cm Radiopaque PEBA

1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm)
(Straight Soft Tip on Both Ends)

1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW

Catheter Priming Volume: 0.16 mL

Priming Volume is approximate.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

⚠ California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

SZ-00522-103B (3/18)

LOT



REF SAC-00522-PBX

22 Ga. | 5 cm catheter length

Arterial Catheterization Set

DO NOT BEND

ARROW