

DO NOT BEND

ARROW



REF SAC-00818-PBX

18 Ga. | 8 cm catheter length

### Arterial Catheterization Set

- cs** Arteriální katetrizační set
- da** Arterielt kateteriserings sæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:  
 1: ..... Indwelling Catheter: 18 Ga. x 8 cm Radiopaque PEBA  
 1: ..... Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)  
 1: ..... Introducer Needle: 18 Ga. x 2" (5 cm) XTW

Catheter Priming Volume: 0.19 mL Priming Volume is approximate.

Rx only.  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65  
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.



EC REP Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
 0086

Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex**  
 www.teleflex.com  
 SZ-00818-105B (3/18)

LOT



DO NOT BEND

ARROW



REF SAC-00818-PBX

18 Ga. | 8 cm catheter length

### Arterial Catheterization Set

- cs** Arteriální katetrizační set
- da** Arterielt kateteriserings sæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:  
 1: ..... Indwelling Catheter: 18 Ga. x 8 cm Radiopaque PEBA  
 1: ..... Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)  
 1: ..... Introducer Needle: 18 Ga. x 2" (5 cm) XTW

Catheter Priming Volume: 0.19 mL Priming Volume is approximate.

Rx only.  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65  
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.



EC REP Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
 0086

Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex**  
 www.teleflex.com  
 SZ-00818-105B (3/18)

LOT



REF SAC-00818-PBX

18 Ga. | 8 cm catheter length

### Arterial Catheterization Set

DO NOT BEND

ARROW



REF SAC-00818-PBX

18 Ga. | 8 cm catheter length

### Arterial Catheterization Set

DO NOT BEND

ARROW