

Dimensions:
Width: 4" (10.2 cm)
Length: 11" (27.9 cm)

DO NOT BEND

REF SAC-00822

Arterial Catheterization Set

ARROW®

- cs** Arteriální katetrizační set
- da** Arterielt kateteriseringssæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:

- 1: Indwelling Catheter: 22 Ga x 3-1/4" (8 cm) Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm)
(Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 22 Ga x 1-1/2" (4 cm) XTW

Catheter Priming Volume: 0.17 ml

Priming Volume is approximate.

Contraindications: Inadequate collateral arterial circulation.

Fluid path components are non-pyrogenic.



Arrow International, Inc., Subsidiary of Teleflex Incorporated
2400 Bernville Road, Reading, PA 19605 USA

EC REP Teleflex Medical, IDA Business and Technology Park
Dublin Road, Athlone, Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

LOT



SZ-00822-101D (3/17)

ARROW®



22 Ga 8 cm catheter length

DO NOT BEND

DO NOT BEND

REF SAC-00822

Arterial Catheterization Set

ARROW®

- cs** Arteriální katetrizační set
- da** Arterielt kateteriseringssæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkatetreringsset
- tr** Arteriyel Kateterizasyon Seti

Contents:

- 1: Indwelling Catheter: 22 Ga x 3-1/4" (8 cm) Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm)
(Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 22 Ga x 1-1/2" (4 cm) XTW

Catheter Priming Volume: 0.17 ml

Priming Volume is approximate.

Contraindications: Inadequate collateral arterial circulation.

Fluid path components are non-pyrogenic.



Arrow International, Inc., Subsidiary of Teleflex Incorporated
2400 Bernville Road, Reading, PA 19605 USA

EC REP Teleflex Medical, IDA Business and Technology Park
Dublin Road, Athlone, Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

LOT



SZ-00822-101D (3/17)

ARROW®



22 Ga 8 cm catheter length

DO NOT BEND