

Dimensions:
Width: 4" (10,2 cm)
Length: 11" (27,9 cm)

DO NOT BEND

REF **SAC-01222**

Arterial Catheterization Set ARROW®

- cs Arteriální katetrizační set
- da Arterielt kateteriseringssæt
- nl Arteriële katheterisatieset
- et Arterikanüüli komplekt
- fi Valtimokatetrisarja
- fr Jeu de cathétérisme artériel
- de Besteck zur Katheterisierung von Arterien
- el Σet αρτηριακού καθετηριασμού
- hu Artériás katéter készlet
- it Set per cateterismo arterioso
- lt Arterijų kateterizavimo rinkinys
- no Arteriell kateteriseringssett
- pl Zestaw do cewnikowania tętnic
- pt Conjunto de cateterização arterial
- ro Set de cateterizare arterială
- ru Набор для артериальной катетеризации
- sk Arteriálna katetrizačná sada
- sl Pribor za arterijsko kateterizacijo
- es Conjunto de cateterización arterial
- sv Artärkateteriseringsset
- tr Arteriyel Kateterizasyon Seti

Contents:

- 1: Indwelling Catheter: 22 Ga x 5" (12 cm) Radiopaque Polyurethane
1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm)
(Straight Soft Tip on Both Ends)
1: Introducer Needle: 22 Ga x 1-1/2" (4 cm) XTW

Catheter Priming Volume: 0.18 ml

Priming Volume is approximate.

Contraindications: Inadequate collateral arterial circulation.

Fluid path components are non-pyrogenic.



Arrow International, Inc., Subsidiary of Teleflex Incorporated
2400 Bernville Road, Reading, PA 19605 USA

EC REP Teleflex Medical, IDA Business and Technology Park
Dublin Road, Athlone, Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

LOT



SZ-01222-101D (3/17)

ARROW®



22 Ga 12 cm
catheter
length

DO NOT BEND

DO NOT BEND

REF **SAC-01222**

Arterial Catheterization Set ARROW®

- cs Arteriální katetrizační set
- da Arterielt kateteriseringssæt
- nl Arteriële katheterisatieset
- et Arterikanüüli komplekt
- fi Valtimokatetrisarja
- fr Jeu de cathétérisme artériel
- de Besteck zur Katheterisierung von Arterien
- el Σet αρτηριακού καθετηριασμού
- hu Artériás katéter készlet
- it Set per cateterismo arterioso
- lt Arterijų kateterizavimo rinkinys
- no Arteriell kateteriseringssett
- pl Zestaw do cewnikowania tętnic
- pt Conjunto de cateterização arterial
- ro Set de cateterizare arterială
- ru Набор для артериальной катетеризации
- sk Arteriálna katetrizačná sada
- sl Pribor za arterijsko kateterizacijo
- es Conjunto de cateterización arterial
- sv Artärkateteriseringsset
- tr Arteriyel Kateterizasyon Seti

Contents:

- 1: Indwelling Catheter: 22 Ga x 5" (12 cm) Radiopaque Polyurethane
1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm)
(Straight Soft Tip on Both Ends)
1: Introducer Needle: 22 Ga x 1-1/2" (4 cm) XTW

Catheter Priming Volume: 0.18 ml

Priming Volume is approximate.

Contraindications: Inadequate collateral arterial circulation.

Fluid path components are non-pyrogenic.



Arrow International, Inc., Subsidiary of Teleflex Incorporated
2400 Bernville Road, Reading, PA 19605 USA

EC REP Teleflex Medical, IDA Business and Technology Park
Dublin Road, Athlone, Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

LOT



SZ-01222-101D (3/17)

ARROW®



22 Ga 12 cm
catheter
length

DO NOT BEND