

DO NOT BEND

ARROW



REF SAC-02318-PBX

18 Ga. | 23 cm catheter length

Arterial Catheterization Set

- cs Arteriální katetrizační set
da Arterielt kateteriseringssett
nl Arteriële katheterisatieset
et Arterikanüüli komplekt
fi Valtimokatetrisarja
fr Jeu de cathétérisme artériel
de Besteck zur Katheterisierung von Arterien
el Σετ αρτηριακού καθετηριασμού
hu Artériás katéter készlet
it Set per cateterismo arterioso
lt Arterijų kateterizavimo rinkinys
no Arteriell kateteriseringssett
pl Zestaw do cewnikowania tętnic
pt Conjunto de cateterização arterial
ro Set de cateterizare arterială
ru Набор для артериальной катетеризации
sk Arteriálna katetrizačná sada
sl Pribor za arterijsko kateterizacijo
es Conjunto de cateterización arterial
sv Artärkatetreringsset
tr Arteriyel Kateterizasyon Seti

Contents:
1: Indwelling Catheter: 18 Ga. x 23 cm Radiopaque PEBA
1: Spring-Wire Guide: .025" (0.64 mm) dia. x 23-5/8" (60 cm)
(Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW

Catheter Priming Volume: 0.28 mL Priming Volume is approximate.

Rx only.
Warning: Read all package insert warnings, precautions, and instructions prior to use.
Failure to do so may result in severe patient injury or death.

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Not made with natural rubber latex.
Fluid path components are non-pyrogenic.



EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
SZ-02318-103B (1/18)

LOT



REF SAC-02318-PBX

18 Ga. | 23 cm catheter length

Arterial Catheterization Set

DO NOT BEND

ARROW

DO NOT BEND

ARROW



REF SAC-02318-PBX

18 Ga. | 23 cm catheter length

Arterial Catheterization Set

- cs Arteriální katetrizační set
da Arterielt kateteriseringssett
nl Arteriële katheterisatieset
et Arterikanüüli komplekt
fi Valtimokatetrisarja
fr Jeu de cathétérisme artériel
de Besteck zur Katheterisierung von Arterien
el Σετ αρτηριακού καθετηριασμού
hu Artériás katéter készlet
it Set per cateterismo arterioso
lt Arterijų kateterizavimo rinkinys
no Arteriell kateteriseringssett
pl Zestaw do cewnikowania tętnic
pt Conjunto de cateterização arterial
ro Set de cateterizare arterială
ru Набор для артериальной катетеризации
sk Arteriálna katetrizačná sada
sl Pribor za arterijsko kateterizacijo
es Conjunto de cateterización arterial
sv Artärkatetreringsset
tr Arteriyel Kateterizasyon Seti

Contents:
1: Indwelling Catheter: 18 Ga. x 23 cm Radiopaque PEBA
1: Spring-Wire Guide: .025" (0.64 mm) dia. x 23-5/8" (60 cm)
(Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW

Catheter Priming Volume: 0.28 mL Priming Volume is approximate.

Rx only.
Warning: Read all package insert warnings, precautions, and instructions prior to use.
Failure to do so may result in severe patient injury or death.

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Not made with natural rubber latex.
Fluid path components are non-pyrogenic.



EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
SZ-02318-103B (1/18)

LOT



REF SAC-02318-PBX

18 Ga. | 23 cm catheter length

Arterial Catheterization Set

DO NOT BEND

ARROW