

**REF AD-04301**

# Central Venous Catheterization Set with Blue FlexTip® Catheter

# ARROW®

**Contents:**

- 1: ..... Indwelling Catheter: 16 Ga x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: ..... Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Luer-Slip Syringe
- 1: ..... Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter

All components are CE 0086 unless otherwise noted.

Rx only.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Keep away from sunlight. Keep dry. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**LOT**

**Use By**

- FR** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip  
Conserver à l'abri du rayonnement solaire. Conserver au sec.
- DE** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter  
Vor Sonnenlicht schützen. Trocken halten.
- IT** Set per cateterismo della vena centrale con catetere Blue FlexTip  
Tenere lontano dalla luce solare. Conservare in un luogo asciutto.
- PT** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip  
Manter afastado da luz solar. Manter seco.
- ES** Conjunto de cateterización venosa central con catéter Blue FlexTip  
Mantener alejado de la luz del sol. Mantener en un lugar seco.
- SV** Set för central venkatetrering med Blue FlexTip kateter  
Skydda från solljus. Förvaras torrt.

**Catheter Diameter:** 16 Ga, 1.7 mm

**Priming Volume:** 0.43 ml  
Priming volumes are approximate and are done without the injection cap.

**Flow Rate:** 2900 ml/hour  
Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

PMS 193

SZ-04301-103F (9/16)

PMS 193