

Dimensions:
 Width: 6-3/4" (17,2 cm)
 Length: 11-5/16" (28,7 cm)

PMS 193

REF ES-04301-P

**Central Venous
 Catheterization Set
 with Blue FlexTip® Catheter**

ARROW®

Contents:

- 1: Indwelling Catheter: 16 Ga x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Pressure Transduction Probe

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use — do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Keep away from sunlight. Keep dry. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
 2400 Bemville Road
 Reading, PA 19605 USA

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

- FR** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
 Conserver à l'abri du rayonnement solaire. Conserver au sec.
- DE** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
 Vor Sonnenlicht schützen. Trocken halten.
- IT** Set per cateterismo della vena centrale con catetere Blue FlexTip
 Tenere lontano dalla luce solare. Conservare in un luogo asciutto.
- PT** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
 Manter afastado da luz solar. Manter seco.
- ES** Conjunto de cateterización venosa central con catéter Blue FlexTip
 Mantener alejado de la luz del sol. Mantener en un lugar seco.
- SV** Set för central venkatetrering med Blue FlexTip kateter
 Skydda från solljus. Förvaras torrt.

Catheter Diameter: 16 Ga, 1.7 mm

Priming Volume: 0.43 ml
 Priming volumes are approximate and are done without the injection cap.

Flow Rate: 2900 ml/hr
 Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



LOT

Use By

5Z04301-107F (9/16)

PMS 193

ARROW
 INTERNATIONAL

REF ES-04301-P 16 Ga 20 cm catheter length .032 inch dia. spring-wire guide
 Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open