

Dimensions:
 Width: 6-3/4" (17,2 cm)
 Length: 11-5/16" (28,7 cm)

PMS 193

REF ES-04730-EK

**Central Venous
 Catheterization Set
 with Blue FlexTip® Catheter**

ARROW®

Contents:

- 1: Indwelling Catheter: 14 Ga x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Luer-Slip Syringe
- 1: Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

Manufacturer
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

EC REP
 CE 0086
 Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
 www.teleflex.com

- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- de** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip
- sv** Set för central venkatetrisering med Blue FlexTip kateter

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.75 ml
 Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

Flow Rate: 5300 ml/hour
 Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- Caution
- Consult instructions for use
- Do not reuse
- Do not resterilize
- Sterilized by ethylene oxide
- Keep away from sunlight
- Keep dry
- Do not use if package is damaged
- Not made with natural rubber latex

LOT

Use By

PMS 193

SL-04730-104G (5/16)

PMS 193

ARROW
 INTERNATIONAL



REF ES-04730-EK **14** Ga **30** cm catheter length **.032** inch dia. spring-wire guide
 Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open