

Dimensions:
Width: 8-5/8" (21.9 cm)
Length: 15-1/2" (39.4 cm)

PMS 259

REF **SI-11142-S**

MAC™
Two-Lumen Central Venous Access Set
with Integral Hemostasis Valve
for use with 7 - 8 Fr Catheters

ARROW®

Contents:

- 1:Two-Lumen Central Venous Access Device:
4" (10 cm) Radiopaque Polyurethane
with 9 Fr Distal Lumen,
Integral Hemostasis Valve,
Extension Line Clamps
- 1:Spring-Wire Guide, Marked: .035"
(0.89 mm) dia. x 17-3/4" (45 cm)
(Straight Soft Tip on One End – "J" Tip on Other)
with Arrow Advancer™
- 1:Introducer Needle: 18 Ga x 2-1/2"
(6.35 cm) XTW
- 1:Arrow® Raulerson Spring-Wire
Introduction Syringe: 5 ml
- 1:Tissue Dilator: 9 Fr (3.0 mm) x 20.7 cm
- 3:Gauze Pad: 4" x 4" (10 x 10 cm)
- 1:Obturator: 8 Fr
- 1:Scalpel: #11
- 1:Cath-Gard® Catheter
Contamination Shield: 31-1/2" (80 cm)
with Tuohy-Borst Adapter, Sterile Tape
- 1:Dust Cap: Vented
- 1:Dust Cap: Non-Vented
- 1:Pressure Transduction Probe
- 1:Stopcock: 4-way High-Flow
- 1:Injection Site Cap

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not sterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Keep away from sunlight. Keep dry. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

LOT

Use
By



(01) 0 08 01902 03290 4

- FR** Jeu d'accès veineux central à deux lumières MAC avec valve hémostatique intégrée, à utiliser avec les cathéters de 7 - 8 Fr
Conserver à l'abri du rayonnement solaire. Conserver au sec.
- DE** Einführbesteck für zweilumigen, MAC-zentralen Venenkatheter mit Hämostaseventil mit Kathetern der Größen 7 - 8 Fr
Vor Sonnenlicht schützen. Trocken halten.
- IT** Set di accesso venoso centrale a due lumi MAC con valvola elastostatica incorporata per l'uso con cateteri da 7 - 8 Fr
Tenere lontano dalla luce solare. Conservare in un luogo asciutto.
- PT** Conjunto para Acesso Venoso Central MAC com Lúmen Duplo, Válvula de Hemostase Integral, para utilizar com Cateters de 7 - 8 Fr
Manter afastado da luz solar. Manter seco.
- ES** Conjunto para acceso venoso central de dos luces MAC con válvula hemostática integral para uso con catéteres de entre 7 y 8 Fr
Mantener alejado de la luz del sol. Mantener en un lugar seco.
- SV** MAC tvåkanaligt set för centralt ventillträde med integrerad hemostasventil för användning med 7 - 8 Fr katetrar
Skydda från solljus. Förvaras torrt.

Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (9 Fr) no catheter	1.65	33,000
Distal (9 Fr) with 8 Fr catheter	0.90	10,500
Proximal (12 Ga) no catheter	0.60	13,000

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

PMS 259

SI-11142-105E (12/14)

PMS 259

ARROW®
INTERNATIONAL

REF **SI-11142-S** 9^{Fr} 2 Lumen 10^{cm} device length .035^{inch dia.} spring-wire guide
MAC™ Two-Lumen Central Venous Access Set

Peel to Open