

Dimensions:  
Width: 6-3/4" (17,2 cm)  
Length: 11-5/16" (28,7 cm)

PMS 259

REF **BR-12703**

Multi-Lumen Central  
Venous Catheterization Set  
with Blue FlexTip® Catheter

**ARROW**®

**Contents:**

- 1: .....Multi-Lumen Indwelling Catheter:  
7 Fr x 6" (16 cm) Radiopaque  
Polyurethane with Blue FlexTip®,  
Extension Line Clamps
- 1: ..... Spring-Wire Guide, Marked:  
.032" (0.81 mm) dia. x 17-3/4"  
(45 cm) (Straight Soft Tip on One  
End – "J" Tip on Other)  
with Arrow Advancer™
- 1: ..... Introducer Needle: 18 Ga  
x 2-1/2" (6.35 cm) XTW  
and 5 ml Luer-Slip Syringe
- 1: ..... Tissue Dilator: 8.5 Fr (2.8 mm)  
x 10.2 cm
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Scalpel: #11
- 3: ..... Stopcock: 4-way High-Flow

- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- PL** Zestaw wielokanałowy do cewnikowania żył centralnych z cewnikiem typu Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetризациjo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkatetrising med Blue FlexTip kateter
- TR** Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile

All components are CE 0086 unless otherwise noted.

Rx only.


Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Keep away from sunlight. Keep dry. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.39	3400
Medial (18 Ga)	0.35	1800
Proximal (18 Ga)	0.37	1900

\*Priming volume are approximate and done without the stopcock.  
†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT

Use By



(01) 0 08 01902 01071 1

**ARROW**  
INTERNATIONAL

REF **BR-12703** 7 Fr 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide  
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open

PMS 259

52-12703-1070 (11/11/14)

PMS 259