

Dimensions:
 Width: 6-3/4" (17,2 cm)
 Length: 11-5/16" (28,7 cm)

PMS 259

REF AD-14702-EK

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1:Two-Lumen Indwelling Catheter:
 7 Fr x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1:Spring-Wire Guide, Marked:
 .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1:Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Luer-Slip Syringe
- 1:Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1:SecondSite™ Adjustable Hub: Fastener
- 1:SecondSite™ Adjustable Hub: Catheter Clamp
- 1:Scalpel: #11
- 2:Stopcock: 4-way High-Flow

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.


Sterile. Fluid path components are non-pyrogenic.

PMS 259

Manufacturer
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

CE REP
 Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
 www.teleflex.com



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.58	3200
Proximal (16 Ga)	0.52	2200

*Priming volume are approximate and done without the stopcock.
 †Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

57-14702-1026 (5/16)

- Caution
- Consult instructions for use
- Do not reuse
- Do not resterilize
- Sterilized by ethylene oxide
- Keep away from sunlight
- Keep dry
- Do not use if package is damaged
- Not made with natural rubber latex


LOT

Use By

PMS 259

PMS 259

ARROW
 INTERNATIONAL



REF AD-14702-EK 7 Fr 2 Lumen 30 cm catheter length .032 inch dia. spring-wire guide
Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open