
Lumen	Priming Vol. <sup>1</sup> (cc)	Flow Rate <sup>2</sup> (cc/hr)
Distal (16 Ga.)	0.58	3200
Proximal (16 Ga.)	0.52	2200

<sup>1</sup>Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 cc.

<sup>2</sup>Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



## Contents:

One: ....Two-Lumen Indwelling Catheter: 7 Fr. x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip® Extension Line Clamps, Injection Site Caps One: . . Spring-Wire Guide, Marked: .032 (.81 mm) dia. x 26-13/16" (68 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™\* ..... Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle One: . . . . . . Fastener: Catheter Clamp One: . . . . . Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW ....Pressure Transduction Probe One: . . . . . Arrow Raulerson Spring-Wire Introduction Syringe<sup>†</sup>: 5 cc . . . . Tissue Dilator One: \* US Patent No. 5,484,419.

## Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

† Licensed under US Patent Nos. 4,813,938 and 5,045,065.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use — do not resterilize. Sterilized by ethylene oxide.


This product contains <u>NO</u> natural rubber latex.

Arrow International, Inc. 2400 Bernville Road Reading, PA 19605 USA

## Product No. ES-14702-J

Two-Lumen
Central Venous
Catheterization Set
with Blue FlexTip® Catheter

## ARROW\*





Product No. **ES-14702-J** 

7<sup>Fr.</sup> 2<sup>Lume</sup>

30 cm catheter length

.032 inch dia. spring-wire guide

