



REF SS-14702

**Two-Lumen CVC for use only with  
Arrow PSI (8.5 Fr. and 9 Fr.)**

**Contents:**


- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) Catheter Contamination Guard
- 2: Dust Cap: Non-Vented

Rx only.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.35	5110
Proximal (16 Ga.)	0.38	4700

\* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an *in vitro* experiment using blood analog solution.

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
www.teleflex.com

SZ-14702-118D (2021-08)

LOT



REF SS-14702

**Two-Lumen CVC for use only with  
Arrow PSI (8.5 Fr. and 9 Fr.)**

**Contents:**


- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) Catheter Contamination Guard
- 2: Dust Cap: Non-Vented

Rx only.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.35	5110
Proximal (16 Ga.)	0.38	4700

\* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an *in vitro* experiment using blood analog solution.

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
www.teleflex.com

SZ-14702-118D (2021-08)

LOT

