



REF **SS-14703**

Three-Lumen CVC for use only with Arrow PSI (8.5 Fr. and 9 Fr.)

Contents:

- 1: Three-Lumen Catheter: 7 Fr. (2.5 mm OD) Contamination Guard (14 cm Catheter Length when used with Percutaneous Sheath Introducer)
- 3: Dust Cap: Non-Vented

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



REF **SS-14703**

Three-Lumen CVC for use only with Arrow PSI (8.5 Fr. and 9 Fr.)

Contents:

- 1: Three-Lumen Catheter: 7 Fr. (2.5 mm OD) Contamination Guard (14 cm Catheter Length when used with Percutaneous Sheath Introducer)
- 3: Dust Cap: Non-Vented

Rx only.


Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.




Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.34	3510
Medial (18 Ga.)	0.27	1690
Proximal (18 Ga.)	0.27	1810

* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an *in vitro* experiment using blood analog solution.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com
SZ-14703-141D (2021-08)



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.34	3510
Medial (18 Ga.)	0.27	1690
Proximal (18 Ga.)	0.27	1810

* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an *in vitro* experiment using blood analog solution.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com
SZ-14703-141D (2021-08)

