

Dimensions:
 Width: 6-9/64" (15,6 cm)
 Length: 8-21/32" (22 cm)

PMS 193

Contents:

- 1: Indwelling Catheter: 16 Ga x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg^{ard}® Antimicrobial Surface Treatment¹, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
 - 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
 - 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
 - 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
 - 1: Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
 - 1: SecondSite™ Adjustable Hub: Fastener
 - 1: SecondSite™ Adjustable Hub: Catheter Clamp
 - 1: Pressure Transduction Probe
- ¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Contraindications: The ARROWg^{ard} Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Sterile. Fluid path components are non-pyrogenic.

Manufacturer
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bemville Road
 Reading, PA 19605 USA

CE REP
 Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
 www.teleflex.com



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex



Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

LOT

Use By

Catheter Diameter: 16 Ga, 1.7 mm

Priming Volume: 0.43 ml

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

Flow Rate: 3650 ml/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip ARROWg^{ard} Blue. Conservez le produit à une température inférieure à 25°C. Veillez à ne pas exposer le produit à une chaleur excessive supérieure à 40°C.
- de** Besteck für zentralen Venenkatheter mit Blue FlexTip ARROWg^{ard} Blue Katheter. Unter 25°C lagern. Zu hohe Temperaturen über 40°C vermeiden.
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip ARROWg^{ard} Blue. Conservare ad una temperatura inferiore a 25°C. Evitare il calore eccessivo superiore a 40°C.
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWg^{ard} Blue. Przechowywać w temperaturze poniżej 25°C. Chronić przed temperaturą powyżej 40°C.
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip, ARROWg^{ard} Blue. Conservar a temperaturas inferiores a 25°C. Evitar temperaturas superiores a 40°C.
- ru** Набор для катетеризации центральных вен с катетером Blue FlexTip ARROWg^{ard} Blue. Хранить при температуре ниже 25°C. Избегать нагрева выше 40°C.
- sl** Komplet za centralno venosko katetrizacijo s katetrom Blue FlexTip ARROWg^{ard} Blue. Shranjujte pod 25°C. Izogibajte se prekomerni vročini nad 40°C.
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip ARROWg^{ard} Blue. Almacenar a temperaturas inferiores a 25°C. Evite temperaturas excesivas por encima de 40°C.
- sv** Set för central venkatetrisering med Blue FlexTip ARROWg^{ard} Blue kateter. Förvaras vid temperatur under 25°C. Undvik höga temperatur över 40°C.
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWg^{ard} Blue Kateter ile 25°C'ın altında saklayın. 40°C'ın üzerindeki aşırı ısıdan sakının.

5Z24301-101H (9/16)

PMS 193

PMS 193

ARROW®



REF CV-24301-E

16^{Ga} 20 cm catheter length .032 inch dia. spring-wire guide

Central Venous Catheterization Set with ARROWg^{ard} Blue® Catheter

Peel to Open