

ARROW®



Arrow-Howes™ Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF **CS-25854-E** | **8.5** Fr. | **4** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

- fr** Jeu de cathétérisme veineux central à quatre lumières Arrow-Howes avec cathéter Blue FlexTip ARROWgard Blue
- de** Arrow-Howes-Besteck für vierlumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWgard-Blue-Katheter
- it** Set per cateterismo della vena centrale a quattro lumi Arrow-Howes con catetere Blue FlexTip ARROWgard Blue
- pt** Conjunto de Cateterização Venosa Central Arrow-Howes com Lúmen Quádruplo e Cateter Blue FlexTip, ARROWgard Blue
- es** Conjunto de cateterización venosa central de cuatro luces Arrow-Howes con catéter Blue FlexTip ARROWgard Blue
- sv** Arrow-Howes set för frykanalig central venkatetrisering med Blue FlexTip ARROWgard Blue-kateter

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.



STERILE EO



Not made with natural rubber latex




Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: The ARROWgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.43	2700
Medial 1 (14 Ga.)	0.54	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.41	1600

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-25854-130C (01/18)

DuPont™ Tyvek®
padding



Arrow-Howes™ Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

ARROW®

REF **CS-25854-E** | **8.5** Fr. | **4** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide