

ARROW®



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF **CS-27702-E** | **7** Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWgard Blue
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWgard-Blue-Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWgard Blue
- pl** Dwukanałowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWgard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWgard Blue
- ru** Комплект для катетеризации главных вен с двухпросветным катетером Blue FlexTip ARROWgard Blue
- sl** Komplet za centralno venosko katetrizacijo z dvema svetlinama s katetrom Blue FlexTip ARROWgard Blue
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip ARROWgard Blue
- sv** Set för tvåkanalig central venkatetrisering med Blue FlexTip ARROWgard Blue-kateter
- tr** İki Lümenli Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWgard Blue Kateter ile



Not made with natural rubber latex

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).


Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: The ARROWgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).




Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.57	5000
Proximal (18 Ga.)	0.40	1500

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

SZ-27702-200A (05/18)



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF **CS-27702-E**

7 Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

ARROW®