

# ARROW®



## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF **CS-27702-E**

**7** Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWgard Blue
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWgard-Blue-Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWgard Blue
- pl** Dwukanałowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWgard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWgard Blue
- ru** Комплект для катетеризации главных вен с двухпросветным катетером Blue FlexTip ARROWgard Blue
- sl** Komplet za centralno venosko katetrizacijo z dvema svetlinama s katetrom Blue FlexTip ARROWgard Blue
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip ARROWgard Blue
- sv** Set för tvåkanalig central venkatetrering med Blue FlexTip ARROWgard Blue-kateter
- tr** İki Lümenli Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWgard Blue Kateter ile



Not made with natural rubber latex


Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Fluid path components are non-pyrogenic.

**Contraindications: The ARROWgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.57	4739
Proximal (18 Ga.)	0.40	880

\* Priming volumes are approximate and are done without accessories.  
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

**Manufacturer:**  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA



SZ-27702-300A (11/19)



## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF **CS-27702-E**

**7** Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

