



Intraosseous Access Tray + Single-Use Power Driver

EZ-IO 15mm Procedure Tray

SKU / Article #: 9018T-VC-005

Considerations

Not made with natural rubber latex.
‡MR Conditional

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Intraosseous Access Tray + Single-Use Power Driver	
SKU / Article #: 9018T-VC-005	
Quantity	Component Description
1	EZ-IO™ Power Driver (Single-Use)
1	EZ-IO™ 15mm 15 Ga Intraosseous Needle‡ Set
1	10 mL Luer-Lock Pre-Filled Saline Syringe
1	1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
1	SharpsAway® II Locking Disposal Cup
1	EZ-IO™ Wristband
1	EZ-Stabilizer™ Dressing
1	EZ-Connect™ Extension Set (Priming Volume - Approx. 1.0 ml)

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



Product Information

ARROW EZ-IO Intraosseous Access Tray

INDICATIONS: For use in the femur, humerus, and tibia of adult patients for the delivery of fluids and medications.

CONTRAINDICATIONS: Do not use in patients with known or suspected bone infection, fracture, or other bone pathology.

WARNINGS: Do not use in children. Do not use in patients with known or suspected bone infection, fracture, or other bone pathology.

PRECAUTIONS: Use aseptic technique. Do not reuse the device.

INSTRUCTIONS FOR USE: See the separate instruction sheet for detailed instructions.

Item	Qty	Part Number	Lot Number	Expiration Date
1	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
2	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
3	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
4	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
5	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
6	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
7	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
8	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
9	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025
10	1	ARROW EZ-IO Intraosseous Access Tray	123456789	12/2025

ARROW EZ-IO Intraosseous Access Tray

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.