



Intraosseous Access Tray + Single-Use Power Driver

EZ-IO 15mm Procedure Tray

SKU / Article #: 9018T-VC-005

Considerations

Not made with natural rubber latex.
‡MR Conditional

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Intraosseous Access Tray + Single-Use Power Driver	
SKU / Article #: 9018T-VC-005	
Quantity	Component Description
1	EZ-IO™ Power Driver (Single-Use)
1	EZ-IO™ 15mm 15 Ga Intraosseous Needle‡ Set
1	10 mL Luer-Lock Pre-Filled Saline Syringe
1	1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® ¹ One-Step Frepp® Solution
1	SharpsAway® II Locking Disposal Cup
1	EZ-IO™ Wristband
1	EZ-Stabilizer™ Dressing
1	EZ-Connect™ Extension Set (Priming Volume - Approx. 1.0 ml)

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



ARROW EZ-IO
INTRAOSSUEOUS ACCESS TRAY

TABLE 1: INTRAOSSUEOUS ACCESS TRAY COMPONENTS

Part Number	Description	Quantity
1000000000	Intraosseous Access Tray	1
1000000001	Single-Use Driver	1
1000000002	Single-Use Power Driver	1
1000000003	Single-Use Intraosseous Needle	1
1000000004	Single-Use Intraosseous Needle	1
1000000005	Single-Use Intraosseous Needle	1
1000000006	Single-Use Intraosseous Needle	1
1000000007	Single-Use Intraosseous Needle	1
1000000008	Single-Use Intraosseous Needle	1
1000000009	Single-Use Intraosseous Needle	1
1000000010	Single-Use Intraosseous Needle	1

© 2024 Teleflex Incorporated. All Rights Reserved.

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.