



Representative photo(s): actual product specifications may vary from that shown in the photo

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

PSI SET: 9 FR

SKU / Article #: AD-09903-S

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters SKU / Article #: AD-09903-S	
Quantity	Component Description
1	Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
1	Dust Cap: Non-Vented
1	Stopcock: 4-way High-Flow

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.