

Representative photo(s): actual product specifications may vary from that shown in the photo

One-Lumen PICC Set

PICC SET: 16 GA X 27 1/2IN

SKU / Article #: AJ-01671



Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

One-Lumen PICC Set SKU / Article #: AJ-01671	
Quantity	Component Description
1	One Lumen Catheter: 16 Ga. (1.80 mm OD) x 70 cm, Side Holes, Blue FlexTip®, Contamination Guard
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Scalpel: #11
1	Stopcock: 3-way

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.