



## Radial Artery Catheterization

### RA CATH KIT: 20 GA X 1-3/4"

SKU / Article #: AK-04020-C

### Considerations

Not made with natural rubber latex.  
Includes Lidocaine


**Sales Quantity/Case: 10**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Radial Artery Catheterization SKU / Article #: AK-04020-C	
Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 27 Ga. x 1/2" (1.27 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 3 mL Luer-Lock
1	Foil Pack 70% Alcohol Prep Pad
1	Skin Protectant Prep
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep <sup>®1</sup> One-Step Frepp <sup>®</sup> Solution
1	Wing Clip: Suture
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
1	Scalpel: #11
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup> A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.