

Radial Artery Catheterization with Integral Needle Protection

RA CATH KIT: 20 GA X 1-3/4"

SKU / Article #: AK-04020-SPC

Considerations

Not made with natural rubber latex.
Includes Lidocaine


Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Radial Artery Catheterization with Integral Needle Protection SKU / Article #: AK-04020-SPC	
Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 3 mL Luer-Lock
1	Skin Protectant Prep
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	1.5 mL Applicator 2% CHG and 70% IPA Chloraprep® ¹ One-Step Frepp® Solution
1	SharpsAway® II Locking Disposal Cup
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.