



Radial Artery Catheterization

RA CATH KIT: 20 GA X 1-3/4"

SKU / Article #: AK-04020

Considerations

Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 10

Radial Artery Catheterization SKU / Article #: AK-04020	
Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 3 mL Luer-Lock
1	Foil Pack 70% Alcohol Prep Pad
1	Skin Protectant Prep
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep® ¹ One-Step Frepp® Solution
1	Wing Clip: Suture
1	Towel: 13" x 18" (33 cm x 46 cm)
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.