

## Small Vein Jugular Puncture

### JUGULAR KIT: 20GA X 5 IN

SKU / Article #: AK-04150-E-S

#### Highlights

---

CE Mark

#### Considerations

---

Not made with natural rubber latex.

**Sales Quantity/Case: 5**

---




NO IMAGE  
AVAILABLE

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Small Vein Jugular Puncture  
SKU / Article #: AK-04150-E-S

Quantity	Component Description
1	One-Lumen Catheter: 20 Ga. (1.2 mm OD) x 12 cm
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
2	Syringe: 3 mL Luer-Lock
1	Tissue Dilator: 4 Fr. (1.35 mm) x 7.6 cm
1	Tissue Dilator: 4 Fr. (1.35 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Scalpel: #11
1	Prep Sponge Swab
3	Gauze Pad: 2" x 2" (5 cm x 5 cm)

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

