

One-Lumen CVC

CVC KIT: 16 GA X 12 IN

SKU / Article #: AK-04400

Considerations

Not made with natural rubber latex.


Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

One-Lumen CVC
SKU / Article #: AK-04400

Quantity	Component Description
1	One-Lumen Catheter: 16 Ga. (1.8 mm OD) x 30 cm
1	Spring-Wire Guide: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
1	Adapter: Arrow-Johans™ RAECG
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
2	Prep Sponge Swab
3	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.