



## One-Lumen CVC

CVC KIT: 14 GA X 8 IN

SKU / Article #: AK-04700

### Considerations

Not made with natural rubber latex.



Sales Quantity/Case: 10

*Representative photo(s): actual product specifications may vary from that shown in the photo*

One-Lumen CVC  
SKU / Article #: AK-04700

Quantity	Component Description
1	One-Lumen Catheter: 14 Ga. (2.2 mm OD) x 20 cm, Side Holes
1	Spring-Wire Guide: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
2	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® Disposal Cup
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Scalpel: #11
2	Prep Sponge Swab
3	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.