



PSI for use with 7 - 7.5 Fr. Catheters

PSI KIT: 8.5 FR

SKU / Article #: AK-09802-A

Considerations

Caution: This product contains dry natural rubber.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo




PSI for use with 7 - 7.5 Fr. Catheters

SKU / Article #: AK-09802-A

Quantity	Component Description
1	Single-Lumen Infusion Catheter, SLIC®: 7 Fr. x 16 cm
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe
1	Syringe: 1 mL Luer-Slip
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ¹ One-Step Solution with Hi-Lite Orange™ Tint
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
1	SharpsAway® Disposal Cup
3	Drape: 26" x 18" (66 cm x 46 cm) with adhesive
1	Scalpel: #11
1	Stopcock: 3-way High-Flow
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Nylon with Straight Needle

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.