



## PSI for use with 7.5 - 8 Fr. Catheters

PSI KIT: 9 FR

SKU / Article #: AK-09903-A

### Considerations

Not made with natural rubber latex.

±Precaution: Contains Phthalate: DEHP

Sales Quantity/Case: 10

*Representative photo(s): actual product specifications may vary from that shown in the photo*



PSI for use with 7.5 - 8 Fr. Catheters

SKU / Article #: AK-09903-A

Quantity	Component Description
1	Arrow-Flex <sup>®</sup> Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow <sup>®</sup> Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep <sup>®1</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint
1	Obturator: 8 Fr.
1	Suture: 3-0 Braided Silk with Straight Needle
1	Cath-Gard <sup>®</sup> Catheter Contamination Shield: 80 cm with TwistLock <sup>™</sup> Adapter
1	Dust Cap: Non-Vented
1	SharpsAway <sup>®</sup> Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Scalpel: #11
1	Stopcock: 4-way High-Flow
7	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Tubing±: 20" (50 cm) Extension

<sup>1</sup> A registered trademark of CareFusion or one of its subsidiaries.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.