

Representative photo(s): actual product specifications may vary from that shown in the photo

Two-Lumen CVC

CVC KIT: 2-LUMEN 7 FR X 16 CM

SKU / Article #: AK-12702-J

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 5

Two-Lumen CVC SKU / Article #: AK-12702-J	
Quantity	Component Description
1	Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
1	Pressure Transduction Probe
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
6	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® Disposal Cup

Two-Lumen CVC SKU / Article #: AK-12702-J	
Quantity	Component Description
1	Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
1	Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Scalpel: #11
2	Stopcock: 4-way High-Flow
2	Prep Sponge Swab
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.