

Small Vein Two-Lumen CVC

CVC KIT: 2-LUMEN 4 FR X 13 CM

SKU / Article #: AK-14412

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5


Representative photo(s): actual product specifications may vary from that shown in the photo

Small Vein Two-Lumen CVC

SKU / Article #: AK-14412

Quantity	Component Description
1	Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
1	Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® Disposal Cup
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.