

Representative photo(s): actual product specifications may vary from that shown in the photo

Three-Lumen CVC

CVC KIT: 3-LUMEN 7 FR X 30 CM

SKU / Article #: AK-14703-SP

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 5




Three-Lumen CVC

SKU / Article #: AK-14703-SP

Quantity	Component Description
1	Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm
1	Spring-Wire Guide: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®1} One-Step Solution with Hi-Lite Orange [™] Tint
3	Dust Cap: Non-Vented
1	SecondSite [™] Adjustable Hub: Fastener
1	SecondSite [™] Adjustable Hub: Catheter Clamp
1	SharpsAway [®] II Locking Disposal Cup
1	SharpsAway [®] Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Safety Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.