

Two-Lumen CVC

CVC KIT: 2-LUMEN 8 FR X 20 CM

SKU / Article #: AK-15802-J

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5



Representative photo(s): actual product specifications may vary from that shown in the photo

Two-Lumen CVC

SKU / Article #: AK-15802-J

Quantity	Component Description
1	Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Introducer Needle: 22 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® Disposal Cup
1	Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Scalpel: #11
2	Stopcock: 4-way
2	Prep Sponge Swab
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.