



## Arrowgard Blue® MAC™ Two-Lumen Central Venous Access for use with 7 - 8 Fr. Catheters

### MULTI-LUMEN/PSI KIT

SKU / Article #: AK-21242-CDC

### Highlights

Antimicrobial: Arrowgard Blue\*  
Sharps Safety

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.  
Includes Lidocaine  
±Precaution: Contains Phthalate: DEHP

### Sales Quantity/Case: 10

#### Arrowgard Blue® MAC™ Two-Lumen Central Venous Access for use with 7 - 8 Fr. Catheters SKU / Article #: AK-21242-CDC

| Quantity | Component Description  |
|----------|--|
| 1        | Two-Lumen Central Venous Access Device: 4-1/2" (11.5 cm) Radiopaque Polyurethane with 9 Fr. Distal Lumen, Arrowgard® Antimicrobial Surface Treatment <sup>1</sup> , Integral Hemostasis Valve, Extension Line Clamps |
| 1        | Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer  |
| 1        | Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle  |
| 1        | Injection Needle: 25 Ga. x 1" (2.54 cm)  |
| 1        | Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)   |
| 1        | Pressure Transduction Probe  |
| 1        | Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe  |
| 1        | Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe   |
| 1        | Syringe: 3 mL Luer-Lock  |


Representative photo(s): actual product specifications may vary from that shown in the photo

|   |   |
|---|---|
| Arrowg+ard Blue MAC Two-Lumen Central Venous Access for use with 7 - 8 Fr. Catheters<br>SKU / Article #: AK-21242-CDC |   |
| Quantity  | Component Description   |
| 1   | Tissue Dilator  |
| 1   | 5 mL 1% Lidocaine HCl Solution and Alcohol Prep   |
| 1   | 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® <sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint |
| 1   | Obturator: 8 Fr.  |
| 1   | Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape                  |
| 1   | Dust Cap: Vented  |
| 2   | Dust Cap: Non-Vented  |
| 1   | SharpsAway® II Locking Disposal Cup   |
| 1   | SharpsAway® Disposal Cup  |
| 1   | Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration   |
| 1   | Drape: 41" x 55" with 4" fenestration   |
| 1   | Needle Holder   |
| 1   | Safety Scalpel: #11   |
| 1   | Stopcock: 4-way High-Flow   |
| 2   | Gauze Pad: 2" x 2" (5 cm x 5 cm)  |
| 5   | Gauze Pad: 4" x 4" (10 cm x 10 cm)  |
| 1   | Suture: 3-0 Silk with Curved Needle   |
| 1   | Suture: 3-0 Braided Silk with Straight Needle   |
| 1   | Tubing±: 20" (50 cm) Extension  |

<sup>1</sup> Licensed under US Patent No. 6,706,024.

<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.