

Small Vein Arrowgard Blue® Two-Lumen CVC

CVC KIT: 2-LUMEN 5 FR X 8 CM

SKU / Article #: AK-25502

Highlights

Antimicrobial: Arrowgard Blue*

Considerations

*Chlorhexidine

Not made with natural rubber latex.

Includes Lidocaine

Sales Quantity/Case: 5

NO IMAGE
AVAILABLE



Representative photo(s): actual product specifications may vary from that shown in the photo

Small Vein Arrowgard Blue® Two-Lumen CVC
 SKU / Article #: AK-25502

Quantity	Component Description
1	Arrowgard Blue® Two-Lumen Catheter: 5 Fr. (1.8 mm OD) x 8 cm
1	Spring-Wire Guide: .021" (0.53 mm) dia. x 17-3/4" (45 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other)
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep® ¹ One-Step Solution with Hi-Lite Orange™ Tint
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® Disposal Cup
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.