

## Pressure Injectable Arrowg+ard Blue Plus® Three-Lumen CVC

PI CVC KIT: 3-L 7 FR X 20 CM AGB+

SKU / Article #: AK-45703-PB1A

### Highlights

Antimicrobial: Arrowg+ard Blue Plus\*  
Sharps Safety

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.  
Includes Lidocaine  
±Precaution: Contains Phthalate: DEHP

**Sales Quantity/Case: 5**

Pressure Injectable Arrowg+ard Blue Plus® Three-Lumen CVC SKU / Article #: AK-45703-PB1A	
Quantity	Component Description
1	Arrowg+ard Blue Plus® <sup>1</sup> Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm, Pressure Injectable
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) RW
1	Injection Needle: SafetyGlide™ <sup>2</sup> 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: SafetyGlide™ <sup>2</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Plus Three-Lumen CVC  
 SKU / Article #: AK-45703-PB1A

Quantity	Component Description
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®3</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint
1	SecondSite <sup>™</sup> Adjustable Hub: Fastener
1	SecondSite <sup>™</sup> Adjustable Hub: Catheter Clamp
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	SharpsAway <sup>®</sup> Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
1	Needle Holder
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Flow Rate Information Card with Injection Log
1	Sterile Procedure Sign
1	Medication Label: 1% Lidocaine
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm <sup>®4</sup> 10 cm x 12 cm
1	Dressing: BIOPATCH <sup>®5</sup>
1	Suture: 3-0 Silk with Curved Needle
1	Tubing±: 20" (50 cm) Extension
1	HemoHopper <sup>®</sup> Fluid Receptacle

<sup>1</sup> Licensed under US Patent No. 7,329,412.






<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>4</sup> A registered trademark of 3M Company.

<sup>5</sup> A registered trademark of Johnson & Johnson Corporation.

# Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.