

NO IMAGE  
AVAILABLE

## Arterial Catheterization Kit

### ARTERIAL CATH KIT: 20GA X 5CM

SKU / Article #: ASK-00520-AHS

#### Highlights

Sharps Safety

#### Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Arterial Catheterization Kit





SKU / Article #: ASK-00520-AHS

Quantity	Component Description
1	Catheter: 20 Ga. x 5 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 5 mL Luer-Slip
1	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	Drape: 30" x 36" (76 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Needle Holder
1	Scissors
3	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>2</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
1	Suture: 3-0 Silk with Curved Needle

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of 3M Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

