

## Arterial Catheterization Kit

### ARTERIAL CATH KIT: 20GA X 5CM

SKU / Article #: ASK-00520-WHO

#### Highlights

Sharps Safety

Sales Quantity/Case: 5

*Representative photo(s): actual product specifications may vary from that shown in the photo*





Arterial Catheterization Kit	
SKU / Article #: ASK-00520-WHO	
Quantity	Component Description
1	Catheter: 20 Ga. x 5 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 5 mL Luer-Slip
1	Skin Protectant Prep
1	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	Drape: 30" x 36" (76 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
3	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>2</sup> 10 cm x 12 cm
1	Dressing: STATLOCK® <sup>3</sup> PICC Plus Catheter Stabilization Device

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of 3M Company.

<sup>3</sup> A registered trademark of C. R. Bard, Inc.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.