



## Arterial Catheterization Kit

### ARTERIAL CATH SET: 18 GA X 6-1/2IN

SKU / Article #: ASK-01618-CHUM

#### Highlights

---

Sharps Safety

#### Considerations

---

Not made with natural rubber latex.

±Precaution: Contains Phthalate: DEHP

**Sales Quantity/Case: 5**

---


*Representative photo(s): actual product specifications may vary from that shown in the photo*

Arterial Catheterization Kit  
SKU / Article #: ASK-01618-CHUM

Quantity	Component Description
1	Catheter: 18 Ga. x 16 cm Radiopaque PEBA
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW
1	Injection Needle: SafetyGlide™ <sup>1</sup> 18 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)
1	Syringe: 10 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
1	Needle Holder
1	Safety Scalpel: #11
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Silk with Curved Needle
1	Tubing±: 24" (61 cm) Pressure Monitoring

<sup>1</sup> A trademark of Becton, Dickinson and Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.